1 Integration, exploration, and analysis of high-dimensional single-cell

2 cytometry data using Spectre

3

Thomas Myles Ashhurst^{*#1-3}, Felix Marsh-Wakefield^{*3,4,5}, Givanna Haryono Putri^{*3,6}, Alanna Gabrielle
Spiteri^{3,7}, Diana Shinko^{1,3}, Mark Norman Read^{3,6,8}, Adrian Lloyd Smith^{1,3}, Nicholas Jonathan Cole
King^{1,2,3,7,9}

7

8 ¹Svdnev Cytometry Core Research Facility, Charles Perkins Centre, Centenary Institute and The University of Sydney, NSW, Australia.²Marie Bashir Institute for Infectious Diseases and Biosecurity, The University 9 of Sydney, NSW, Australia. ³Charles Perkins Centre, The University of Sydney, NSW, Australia. ⁴School of 10 11 *Medical Sciences, Faculty of Medicine and Health, The University of Sydney, NSW, Australia.* ⁵*Vascular* Immunology Unit, Department of Pathology, The University of Sydney, NSW, Australia. ⁶School of 12 Computer Science, The University of Sydney, NSW, Australia. ⁷Viral Immunopathology Laboratory, 13 Discipline of Pathology, School of Medical Sciences, Faculty of Medicine and Health.⁸The Westmead 14 *Initiative, Faculty of Engineering, The University of Sydney, NSW, Australia.* ⁹*Sydney Nano, The University* 15 of Sydney, NSW, Australia 16

17

18 *Authors contributed equally to this work

19

20 #Corresponding author:

21 Thomas M. Ashhurst

22 Charles Perkins Centre, The University of Sydney

- 23 Camperdown, NSW, 2050
- **24** + 61 2 8627 1741
- 25 thomas.ashhurst@sydney.edu.au

26

Keywords: high-dimensional cytometry, spectral cytometry, mass cytometry, computational analysis,
clustering, dimensionality reduction, FlowSOM, t-SNE, UMAP

29

30 ABSTRACT

31 As the size and complexity of high-dimensional cytometry data continue to expand, comprehensive, 32 scalable, and methodical computational analysis approaches are essential. Yet, contemporary clustering and 33 dimensionality reduction tools alone are insufficient to analyze or reproduce analyses across large numbers 34 of samples, batches, or experiments. Moreover, approaches that allow for the integration of data across 35 batches or experiments are not well incorporated into computational toolkits to allow for streamlined 36 workflows. Here we present Spectre, an R package that enables comprehensive end-to-end integration and 37 analysis of high-dimensional cytometry data from different batches or experiments. Spectre streamlines the 38 analytical stages of raw data pre-processing, batch alignment, data integration, clustering, dimensionality 39 reduction, visualization and population labelling, as well as quantitative and statistical analysis. Critically, 40 the fundamental data structures used within Spectre, along with the implementation of machine learning 41 classifiers, allow for the scalable analysis of very large high-dimensional datasets, generated by flow 42 cytometry, mass cytometry (CyTOF), or spectral cytometry. Using open and flexible data structures, 43 Spectre can also be used to analyze data generated by single-cell RNA sequencing (scRNAseq) or high-44 dimensional imaging technologies, such as Imaging Mass Cytometry (IMC). The simple, clear, and modular 45 design of analysis workflows allow these tools to be used by bioinformaticians and laboratory scientists

46 alike. Spectre is available as an R package or Docker container. R code is available on Github
47 (<u>https://github.com/immunedynamics/spectre</u>).

48

49 INTRODUCTION

50 High-dimensional analysis tools

51 High-dimensional (HD) cytometry plays an important role in the study of immunology, infectious diseases, 52 autoimmunity, hematology and cancer biology, elucidating critical mediators of immunity and disease at a single-cell level. Advances in single-cell cytometry systems (including flow, spectral and mass cytometry) 53 54 have enabled the simultaneous analysis of over 40 proteins [1] in a single panel, resulting in vast and complex datasets. A large portion of the analysis still utilizes manual gating, which is the sequential and 55 56 often arduous process of identifying cells, based on the expression of one or two cellular markers at a time. While this allows for user-guided examination of the dataset, it is intractable when mapping out the vast 57 variety of cell phenotypes that may be present in HD datasets, and may introduce selective bias through a 58 59 subjective and sequential bifurcating inclusion/exclusion of markers [2]. This is in part due to the number of gates required to fully parse the dataset. As such, a variety of computational approaches have been 60 adopted by the cytometry community to help analyze HD datasets, including automated gating [3], 61 62 clustering (such as PhenoGraph [4], FlowSOM [5], X-Shift [6]), dimensionality reduction (such as t-SNE 63 [7, 8], UMAP [9], trajectory inference (such as Wanderlust [1], Wishbone [10]), and automated cell 64 classification [11] [12, 13]. Many of these tools have been brought together into 'toolboxes', providing either code- or GUI-based analysis workflows, such as Cytofkit [14], CITRUS [15], CATALYST [16], 65 66 Cytofworkflow [17], or diffcyt [18]. In parallel, similar tools have also been developed in other fields, such 67 as a wide variety of analysis approaches for single-cell RNA sequencing (scRNAseq; Seurat [19, 20], 68 SingleCellExperiment [21], Monocle [22]). These have been developed to address similar difficulties in analysis, with developments in one field sometimes assisting another, such as the graph-based clustering 69

from Phenograph [4] informing the design of clustering in Seurat [19]. Overall, these tools allow for
automated and data-driven data processing that may be performed in an unsupervised (requiring no human
supervision) or semi-supervised manner, requiring some human decision-making.

73 Limitations of existing algorithms and toolkits

74 Despite the advantages offered by these computationally-driven approaches, a number of challenges persist, 75 including slow operation speeds, difficulty in scaling to large datasets, and insufficient reproducibility of 76 analysis results across independent experiments. As such, the use of such tools is often limited to relatively small datasets from single experiments. In particular, many scRNAseq specific tools, while suited to 77 processing datasets with a large number of features (markers), do not necessarily scale well to dataset 78 79 volumes that include tens to hundreds of millions of cells, as the volume of data may exceed the memory capacity of the computer being used for analysis, and/or the processing may carry excessively long run 80 81 times [23]. In scRNAseq, independent datasets of this volume are not common outside of cell atlas-type projects [24], but cytometry frequently generates datasets of this size. In addition, these data-driven 82 83 approaches are highly sensitive to the influence of batch effects, a phenomenon in which cells of the same phenotype differ in their signal intensity across multiple experimental batches. Batch effects are an artefact 84 85 of time and context, in which experimental conditions and/or instrumental performance inadvertently 86 influence the measured signal intensity. Left uncorrected, computational algorithms may falsely identify 87 differences (or similarities) based on batch, rather than biologically relevant (and correct) differences between samples or experimental groups. Hence, having the ability to control for and differentiate batch 88 89 effects from biological differences is critical. While a variety of approaches exist for scRNAseq data [19, 20, 25-27], approaches designed for cytometry data are less well developed. While there exist several 90 91 approaches to align data from different batches, many of these align samples individually against each other 92 in the context of standardized clinical profiling [28], potentially removing important biological variations 93 between samples. In contrast, recent approaches calculate batch adjustments using reference controls: 94 aliquots of a control sample run with each batch. The resultant conversions can then be applied to all

95 samples in the batch, preserving biologically relevant differences [29]. Because some batch effects can 96 differentially impact select cell lineages [30], further techniques extend this reference-based approach to 97 apply adjustments in a cluster-specific manner [30]. As this is a developing area, a number of other 98 approaches have also been proposed, using a variety of methods to remove technical variance [11, 31], 99 some inspired by approaches in scRNAseq.

100 While a number of clustering, dimensionality reduction, and batch alignment approaches exist, many of 101 these are not directly integrated with each other - and many implementations require specific dataset 102 formats, limiting the flexibility and interoperability of these tools. Specifically, tools that are developed in 103 different fields (e.g. cytometry or scRNAseq) are often provided as stand-alone packages that operate on 104 custom data formats (e.g. flowFrames, SingleCellExperiment objects, Seurat objects, etc). The dependency 105 of individual tools on custom data formats also makes it non-trivial to apply computational tools from a 106 broader data science and machine learning field that are not specifically designed for cytometry (and thus 107 have no knowledge of the custom data format) but are of significant advantage in existing analysis pipelines. 108 A streamlined solution that can operate on large datasets at speed and integrate tools across fields in a 109 modular fashion would be of significant advantage to the cytometry community.

110 Spectre for analysis of large and complex high-dimensional cytometry datasets

To address these challenges, we developed Spectre, a computational toolkit in R that facilitates rapid and 111 112 flexible analysis of large and complex cytometry datasets. This toolkit expands on the 'Cytometry Analysis Pipeline for large and compleX datasets' (CAPX) workflow that we have published previously for deep 113 profiling of hematopoietic datasets [32]. Through specific function and workflow design, we demonstrate 114 intuitive, modular, and high-speed workflows for data pre-processing (such as scaling/transformation), 115 116 clustering, dimensionality reduction, plotting, as well as quantitative statistical analyses. Additionally, we incorporate an effective method for integrating data across multiple batches or experiments by extending 117 the functionality of CytoNorm [30] batch alignment algorithm. Furthermore, we provide a means to transfer 118 119 labels from reference datasets onto new datasets using machine learning classification algorithms, allowing

for reproducible analytical workflows across multiple experiments. An extension of this process allows for the analysis of very large datasets, where the size of the dataset may exceed the memory capacity of the computer being used. Analysis with Spectre can be applied to flow, spectral, and mass cytometry datasets, consisting of tens to hundreds of millions of cells, as well as to scRNAseq data and HD imaging data (such as Imaging Mass Cytometry, IMC). Spectre is available as an R package or Docker container from Github (https://github.com/immunedynamics/spectre).

Here we demonstrate the utility of the package by analyzing cells isolated from murine bone marrow (BM),
spleen, and central nervous system (CNS) following inoculation with West Nile virus (WNV), measured
by HD flow, spectral, or mass cytometry (CyTOF). Because of the significant cellular dynamics that occur
in multiple tissues in response to infection with WNV, this provides an ideal model for demonstrating
analysis workflows using Spectre.

131

132 METHODS

133 Sample preparation and acquisition

134 Ethical approval for the experimental use of mice was obtained from the Animal Ethics Committee at the 135 University of Sydney. Briefly, mice were anesthetized with 250-300 μ L of Avertin anesthetic via an 136 intraperitoneal injection, and then inoculated intranasally (i.n.) with 10 µL of PBS or a lethal dose of West Nile virus (WNV) in sterile PBS ($6x10^4$ plaque-forming units [PFU]). After 1-7 days post infection (*p.i.*), 137 138 mice received 350 µL of Avertin anesthetic, followed by vena cava section and left ventricular cardiac perfusion with 10 - 30 mL of ice-cold PBS. Spleen, BM and CNS tissue were then extracted and prepared 139 for flow/spectral [33] or mass cytometry [32] as previously described. Labelled samples were acquired on 140 a 5-laser BD LSR-II, a 10-laser custom BD LSR-II, a 5-laser Cytek Aurora, or a Fluidigm Helios. Initial 141 142 cleanup gating and compensation was performed using FlowJo v10.7, where single live leukocytes were 143 then exported as CSV (scale value) or FCS files.

144 Data management in R via data.table

Data management and operations within Spectre were performed using the *data.table* format, and extension of R's base *data.frame*, provided by the *data.table* package [34]. This table-like structure organizes cells (rows) against cellular features or metadata (columns). This simple *data.table* structure allows for the highspeed processing, manipulation (subsetting, filtering, etc.), and plotting of large datasets, as well as fast reading/writing of large CSV files. While CSV files are the preferred file format in this context, Spectre also supports reading and writing FCS files through the use of *flowCore* [35], where the data are then converted from an FCS file into a *flowFrame*, and then into a *data.table*.

152 Data pre-processing and transformation

In Spectre, the entire dataset to be analyzed is condensed into a single *data.table*. Each row (cell) contains cellular expression data and metadata pertaining to the file, sample, group, and batch, stored in separate columns. Metadata for each of the columns (e.g., alternative cellular marker names, etc.,) in the *data.table* can be stored separately, if required.

In order to reduce the contribution of background to measured signal, and to convert cytometry data into a 157 158 linear space, we facilitated ArcSinh transformation of data using the *do.asinh* function in Spectre. This 159 function allows a transformation co-factor to be specified for all columns (typical in mass cytometry) or for individual columns (typical in flow or spectral cytometry). The resulting transformed values are added to 160 161 the dataset as new columns. This addition of new values, rather than replacement of original values, enables 162 greater data retention and transparency. Additional pre-processing steps, such as noise reduction for values below zero, and normalization of values onto the range of zero to one, are provided in the *do.noise.reduce* 163 and *do.normalise* functions respectively, which similarly adds columns of new data. 164

165 Batch alignment and integration

To facilitate batch alignment in Spectre, we developed a wrapper for CytoNorm [<u>30</u>] using the *do.align*function. CytoNorm utilizes reference samples (aliquots of a control sample run with each batch) to

168 determine batch-derived technical differences and calculates a quantile conversion model to align data from 169 each of these batches. This conversion model is then applied to samples in each batch, removing technical 170 variation while preserving biologically relevant differences. The *do.align* function operates on a *data.table*, 171 and allows for reference and target data to be easily specified in a dataset containing a mixture of batches. 172 Quantile conversions for the entire dataset, or conversions for individual clusters, can then be calculated, 173 resulting in a new set of 'aligned' features that are compatible across the batches. Since raw, ArcSinh 174 transformed, and aligned data are preserved in the *data.table*, transparency of these analytical processes is maintained throughout the analytical process. 175

176 Clustering of cells/events

177 Clustering was implemented using the *run.flowsom* function, which acts as a wrapper around the FlowSOM function (available as 'FlowSOM' from Bioconductor [36]), which organizes acquired events (i.e., cells) 178 179 into clusters using a self-organizing map (SOM), and then groups these clusters together into 'metaclusters'. The function accepts a *data.table* and set of column names for generating clusters as input. Once clustering 180 181 is performed, *run.flowsom* returns the *data.table* with new columns containing the cluster and metacluster IDs for each cell. Spectre's run.flowsom provides two options for specifying a target number of 182 183 metaclusters. By default, FlowSOM will determine the number of metaclusters to generate by performing 184 consensus hierarchical clustering on the clusters using a range of different metacluster numbers. It then 185 computes the variance in each metacluster number and determines the point in which the variances suddenly decrease at a gentler rate (the elbow point), which is regarded as the optimal number of metaclusters for the 186 dataset. We refer the reader to the FlowSOM publication for more details on how the optimum number is 187 chosen [5]. Alternatively, users can manually specify a target number of metaclusters (which is 188 189 recommended over automatic determination) or chose not to generate metaclusters at all. In addition to the 190 ability to choose the number of metaclusters to create, the user is also able to change the SOM grid size 191 (xdim and ydim parameters), the seed used to generate clusters and metaclusters for reproducibility.

192 Dimensionality reduction

193 Non-linear dimensionality reduction (DR) was implemented in Spectre using the *run.tsne*, or *run.umap* 194 functions, which are wrappers around t-SNE (available as 'rtsne' from CRAN [37]), and UMAP (available 195 as umap from CRAN [38]), respectively. Both the *run.tsne* and *run.umap* functions accept a *data.table* and set of column names for DR as input, and return the *data.table* with new columns containing tSNE or 196 197 UMAP coordinates for each cell. We also provide a function for running Principal Component Analysis 198 (PCA), a linear DR approach, using the *run.pca* function, acting as a wrapper around *prcomp* (available in the 'stats' package from CRAN [39]). The input data for *run.pca* can either be individual cells (to find cell 199 200 markers that contribute to the variance), or individual samples/patients with summary data (such as median 201 fluorescence intensity (MFI), cell counts or proportions). The output of run.pca is similar to run.tsne and 202 run.umap, though additional outputs can be generated, including scree and contribution plots for detailed 203 assessment of the source of variance in the dataset.

204 Plotting and visualization

205 Plotting of cellular data was implemented in the *make.colour.plot* function, serving as a wrapper around 206 ggplot2 [40]. The dataset, and desired columns to use as X and Y axis are taken as input. By default, each plotted cell is colored by relative plot density. Where desired, each cell can be colored by the level of 207 expression of a specified marker (column), or color by some factor (e.g. cluster, sample, or group etc). To 208 209 reduce the effect of outlier datapoints on the color scale, minimum and maximum thresholds are applied by default as 0.01 (1st percentile) and 0.995 (99.5th percentile) respectively, so that datapoints with values 210 below or above these thresholds are colored at the minimum or maximum respectively. These can be 211 212 modified where required. In cases where a subset of the dataset is being plotted (e.g. separate samples or groups, etc) the limits of the X, Y, and color parameters can be set by the complete dataset, allowing for 213 214 the plots to be compared directly. By default, plots will automatically be saved to the current working 215 directory, allowing for the rapid generation of plots. This plotting functionality is extended in the make.multi.plot function, which arranges multiple plots consisting of set of colored plots (e.g. for examining 216 marker expression) or set of sample group divisions (e.g. for comparing changes between groups). 217

218 Aggregation and summary data

A number of analysis steps require the aggregation or summary of data, either by sample, by cluster, or by 219 220 both. To facilitate this, we implemented two functions - do.aggregate and write.sumtables. The 221 *do.aggregate* function takes advantage of the fast aggregation function of *data.table*, providing the mean 222 or median level of expression of each marker per cluster/population. This is helpful for plotting expression 223 heatmaps, allowing for the interpretation of cluster/population identities. For every population in each 224 sample, the *write.sumtables* function will calculate the percentage of each population as a percentage of cells in each sample (as well as cell counts for each population per sample if a total count of cells per sample 225 is provided). Additionally, for all specified cellular markers, the median expression level, as well as the 226 227 percentage of each population that expresses the marker are calculated for each population in every sample. 228 These summarized *data.tables* are saved to disk as CSV files, in the format of samples (rows) vs features (columns), which can then be used for quantitative and statistical analysis, either in the same workflow, or 229 230 in a separate pipeline.

231 Heatmaps

To facilitate the generation of heatmaps for examining marker expression on populations, or for comparing populations across samples, we implemented '*make.pheatmap*' as a wrapper around the pheatmap function from pheatmap [41]. In *make.pheatmap*, the user must specify the dataset, the column that contains values to be used as the heatmap rows (e.g. sample or cluster identity), and the columns to be plotted as heatmap columns. This can be used to examine expression of marker (columns) on each cluster (rows), or can be used on z-score transformed data plotting immune features (columns) vs samples (rows), to reveal data variance for each cellular feature being measured.

239 Group and volcano plots

To facilitate quantitative and statistical analysis, we implemented functions to create graphs
(*make.autograph*) and volcano plots (*make.volcano.plot*). Graphs are constructed using the *ggplot2* R

242 package, and integrate pairwise statistical comparisons between groups using functionality from the ggpubr 243 R package [42]. Pairwise comparisons between pairs of groups are performed using either Wilcox test (equivalent to a Mann-Whitney test), or a t-test. Overall group variance is assessed using a Kruskal Wallace 244 245 test, or a one-way ANOVA. Volcano plots were generated using the EnchancedVolcano package from 246 Bioconductor [43], and customized to work seamlessly with *data.table* data structures. Initially the relative -fold change of values in each immune parameter between two groups are calculated, along with the p-247 248 value of that change. Each immune parameter can then be plotted as fold change vs p-value, allowing for a rapid assessment of changes with strong and significant effects between groups. 249

250 Classification

251 Following batch alignment, the transfer of population labels from one dataset to another was implemented via two functions: run.train.knn.classifier and run.knn.classifier. The run.train.knn.classifier function uses 252 253 the caret [44] package to train a k-nearest neighbor (kNN) classifier [42] using the labelled dataset, allowing 254 for the identification of the two most closely-matched cells across each dataset. It first normalizes each 255 feature (marker) of the dataset to a range between 0 and 1 using the preprocess function of caret. This is to 256 ensure a uniform distribution of values across all features. Thereafter it trains the KNN classifier using the 257 n-fold cross-validation (n-fold CV) technique. The n-fold CV technique splits the data into n 258 complementary subsets and uses n-1 of those to train the classifier (training data), and the remaining one to 259 test the effectiveness of the trained classifier (testing data). This process is repeated n times, whereby in each iteration different subsets of data are used for training and testing data. For n-fold CV, the overall 260 261 quality of the classifier is computed as the average quality obtained for all n rounds of training and testing. The run.knn.classifier function then runs a kNN classifier provided by FNN R package [45] on a data.table 262 263 without cellular population labels, returning the *data.table* with a new column of population labels added.

264


265 **RESULTS**


266 Spectre facilitates comprehensive end-to-end integration and analysis of large high-dimensional 267 cytometry datasets

Spectre was developed to facilitate rapid and flexible analysis of large and complex cytometry datasets 268 across multiple batches or experiments. Specifically, Spectre facilitates data pre-processing (Figure 1A), 269 270 alignment of data from multiple batches/experiments (Figure 1B), clustering (Figure 1C), dimensionality 271 reduction and visualization (Figure 1C), manual or automated population classification and labelling 272 (Figure 1C), as well as extensive plotting and graphing options for qualitative and quantitative statistical analyses (Figure 1D). Key to this process is strategic selection, implementation, and customization of high-273 274 performance computational tools; and the development of wrapper functions around these tools, enabling 275 them to operate on and produce the same data format for input and output, respectively. This allows for 276 multiple analysis and plotting/graphing tools to be seamlessly weaved together into single analysis workflows, where functions can be used throughout the stages of the workflow, drastically increasing ease 277 278 of use. Moreover, this modularity and flexibility allows these workflows to be adapted to meet different 279 experimental requirements or analytical approaches. Spectre can be used on data generated by a variety of single-cell technologies, including flow, spectral, and mass cytometry (Figure 1E). In addition, following 280 281 some additional pre-processing, Spectre can also be used for the analysis of data generated by scRNAseq 282 or HD imaging technologies such as IMC. (Figure 1E).

283 Data management and pre-processing

Critical to our approach is the use of the *data.table* data structure and operations, in place of prevalent FCS/flowFrame objects or other data structures. The *data.table* format is an enhancement of base R's *data.frame*, a table-like structure commonly used to store data. This enhancement allows for fast and efficient aggregation and manipulation of data through the use of concise flexible syntax and low-level parallelism [34] (Supp. Figure 1). Using Spectre, all samples in the analysis are merged into a single *data.table*, with relevant sample, group, and batch information stored in separate columns. As a result, each








Figure 1. Spectre analysis overview

An overview of Spectre's analysis workflow. A) Data preparation steps including sample, group, and batch annotation, in addition to ArcSinh transformation. B) Batch alignment using CytoNorm. C) Clustering and dimensionality reduction, along with marker expression plotting and expression heatmaps. D) Quantification and statistical analysis through z-score/fold-change heatmaps and grouped dot plots. E) Application of Spectre's analysis workflow to data generated by different technologies, including analyzing a split murine spleen sample by flow cytometry, spectral cytometry, and mass cytometry. Also shown are unrelated PBMC data analyzed by single-cell RNA sequencing (scRNAseq) acquired from 10X genomics, via the Seurat webpage (https://satijalab.org/seurat/v3.2/pbmc3k_tutorial.html) and imaging data generated by Imaging Mass Cytometry (IMC).

row contains all the information relevant for a particular cell, making data manipulation and filtering with
 data.table simple and fast.

A key initial step in computational analysis is the transformation of cellular expression data. Biologically 292 293 meaningful results are most easily interpreted through plotting cellular expression on a logarithmic scale. 294 Because of potentially misleading visual artifacts for signals at the low end of the logarithmic scale, the 295 logicle/bi-exponential scale was developed, where the high end of the scale is logarithmic and the low end 296 of the scale is converted into a linear scale, and the scale then returns to logarithmic at values below the linear component (Supplementary Figure 2A) [46, 47]. Critically, this allows for the compression of low-297 298 end data points with high spreading error, autofluorescence, or noise into a linear space around zero, which 299 can be tailored for the requirements of each channel (Supplementary Figure 2B). For computational analysis 300 to meaningfully manage biological data, a similar compression of low-end data needs to be performed. In 301 is cytometry, this commonly performed using the ArcSinh 302 (http://mathworld.wolfram.com/InverseHyperbolicSine.html) transformation [48]. The data values are 303 transformed into a linearized format, where compression of low-end values is determined using a specified co-factor to determine the extent of compression around zero (Supplementary Figure 2B). 304


Using Spectre, ArcSinh transformation was applied to data from each channel/marker individually, using
different co-factors, allowing for highly customizable data transformations. For flow, spectral, and mass
cytometry data, we tested various co-factor values across multiple channels (Supplementary Figure 3-4).
Overall, we found that a co-factor of 5 to 15, was suitable for all mass cytometry channels. However, in our
experience, we found a range between 100 to 10,000 to be suitable for different channels in conventional
or spectral flow cytometry data.

Because low-end data compression through ArcSinh transformations may still result in significant spread of the negative population around zero, we reasoned that this could be reduced by converting any values that occur below zero, to zero, effectively reducing the contribution of noise below zero. This was implemented in the *do.noise.reduce* function, where a minimum threshold can be set, where all values below that value will be converted to that value (Supplementary Figure 5). Additionally, we developed a function for re-distributing (aka 'normalizing') ArcSinh data between two new values, usually 0 and 1. This prevents markers with extremely high expression levels from exerting greater influence over clustering and DR results, when compared to other markers. This was implemented in the *do.normalise* function (Supplementary Figure 5). Although we have demonstrated the utility of the *do.noise.reduce* and *do.normalise* functions in Supplementary Figure 5, they were not used for the rest of the data presented in this paper.


322 Integrating data from multiple batches or experiments into a single feature space


323 When samples are prepared, stained or run in multiple batches, technical batch-effects can occur, usually 324 consisting of shifts in signal intensity in one or more markers. Because of this, data clustered together with uncorrected batch effects may separate samples based on the batch they belong to, a confounding factor 325 326 that substantially hinders aggregated analysis of datasets from multiple batches or experiments. To provide a comprehensive and adaptable batch alignment and data integration approach, we expanded on the 327 328 functionality of CytoNorm within Spectre. Users specify reference control samples that Spectre uses to determine the alignment conversions. Typically, these are aliquots of a single 'healthy' patient sample that 329 330 are run with each batch of samples. These reference samples should span the full range of the data seeking 331 to be aligned. Where marker expression data is absent on the reference controls (e.g. absence of activation 332 markers on cells from healthy donors), alignment is not feasible, but the original expression data for those markers can still be analyzed. In cases where multiple aliquots of a single control sample are not available, 333 334 multiple control samples of the same type (e.g., healthy mouse BM) that are run with each batch can be used as the next best. Finally, it is also possible to use all samples from each batch as 'reference' controls, 335 336 taking the entire range of marker expression from each batch into account. However, these later approaches should be used with caution, as any biological differences between the samples used for reference will be 337 interpreted as technical variation due to batch effects. To execute this alignment process in Spectre, user-338 indicated reference samples are extracted from the combined *data.table* (Figure 2A) and clustered using 339






Minor batch effects





Major batch effects






Figure 2. Batch alignment using CytoNorm

Batch alignment process using CytoNorm. A) Reference samples acquired with each batch are extracted from the *data.table*, and B) clustered using FlowSOM. Metacluster-specific quantile conversion models are then calculated. C) Cells from all samples/batches are mapped to the FlowSOM grid, and assigned to their nearest metacluster. Cells are then aligned using the metacluster-specific quantile conversion models D) Two sets of BM samples with synthetic batch effects introduced in a population-specific manner. E) CytoNorm alignment was performed using a metacluster-specific alignment process (fine alignment). F) Two sets of BM samples generated with slightly different panels, but targeting the same cellular markers, resulting in significant batch effects. G) CytoNorm is initially performed on the whole dataset (coarse alignment) by mapping the entire dataset into a single metacluster, where H) subsequent FlowSOM clustering allowed for further metacluster-specific alignment (fine alignment).

FlowSOM. For each resulting metacluster, quantile conversion coefficients between each batch for each marker are calculated (Figure 2B). Cells from the full dataset are then mapped to the FlowSOM grid, and each cell is assigned to its nearest metacluster (Figure 2C). Quantile conversion is then applied to the cells in each metacluster, unifying cells from each batch into a single feature space, while preserving biological differences between experimental groups (Figure 2C).

345 To verify the robustness of this process, we applied batch alignment to a set of mock- or WNV-infected 346 BM samples, where synthetic data manipulations were introduced in a population-specific manner to half the samples, to mimic population-specific batch effects. One healthy BM sample from each 'batch' were 347 selected as a reference controls, aggregated together, and plotted using UMAP, revealing subsets of 348 349 neutrophils, eosinophils, monocytes, B cells, T/NK cells and progenitors. Batch-specific differences in their 350 distributions were evident in the UMAP plots (Figure 2D). When FlowSOM was run, consistent populations from each batch were captured within the same metacluster, despite the presence of batch effects. As such, 351 352 metacluster-specific alignment with CytoNorm was able to adequately integrate the cells from each 353 population into a unified dataset (Figure 2E).

354 To test this process on datasets with more substantial batch-effects, we applied batch alignment to two sets of mock- or WNV-infected BM samples, prepared and acquired months apart, using slightly different 355 356 panels. After matching column names between the two datasets, the disparity of these two datasets was 357 apparent upon plotting with UMAP (Figure 2F). In this scenario, the batch effects were so significant that cells from the same population in each batch did not map to the same metacluster, thus preventing 358 359 alignment. To address this, we initially performed a 'coarse' alignment by mapping all cells from each 360 batch into a single metacluster, and performing quantile alignment on the entire dataset. This corrected the 361 majority of batch effects, though some metacluster-specific effects were still evident (Figure 2G, inset). Nevertheless, the resulting data was sufficiently aligned so that populations from each batch could now be 362 mapped into consistent metaclusters (Figure 2G), allowing for a more fine-tuned alignment (Figure 2H) 363

that corrected residual metacluster-specific batch effects (Figure 2H, inset) and unified the dataset into asingle feature space.

366 Clustering and dimensionality reduction strategies

367 A critical analytical step in HD cytometry is the comprehensive identification of cellular populations, 368 including those that are well-established and those that are yet to be characterized. This is particularly 369 relevant when comparing between experimental groups (e.g., diseased patients compared to healthy 370 controls). Clustering and dimensionality reduction are powerful techniques for identifying cell populations 371 and contrasting them between groups. Clustering tools collect phenotypically similar cells into groups (clusters) in a data-driven fashion (i.e., cells are grouped together based on the similarity of their marker 372 373 expression). The output of many clustering approaches are plotted as a collection of nodes, where each node represents a cluster that contains a number of cells. These nodes are typically connected by a form of 374 375 minimum spanning tree (MST) [5] and colored by mean or median marker expression of the cells within each cluster. However, verifying the phenotypic heterogeneity (or homogeneity) of cells captured within a 376 377 cluster is difficult when looking at the data at the cluster level. An alternative approach is to compress 378 cellular data onto two dimensions using non-linear dimensionality reduction tools such as t-SNE [7, 8] and 379 UMAP [9], and to visualize these on a scatter plot, coloring each cell by marker expression level or cluster/metacluster ID. 380

381 Spectre supports clustering and DR by providing a wrapper function for FlowSOM, t-SNE, and UMAP. These functions accept and return data in the *data.table* format. Whilst clustering tools such as FlowSOM 382 scale well to large datasets [49], DR approaches such as t-SNE and UMAP do not; as they incur lengthy 383 computing time, excessive memory usage, and significant crowding effects that inhibit their utility 384 385 (Supplementary Figure 6). Whilst some improvements to runtime (flt-SNE [50]) and plot crowding (opt-386 SNE [51]) have been made, scalability and plot crowding limitations persist. As DR tools are primarily used to visualize cellular data and clustering results, we propose plotting only a subset of the clustered data, 387 which addresses scalability and retains legibility. By using proportional subsampling from each sample, the 388


relative number of cells from each cluster in each sample can be preserved in a smaller dataset, allowing for interpretable analysis via DR. Putative cellular populations amongst the clusters can then be identified, and annotated in both the subsampled DR dataset, as well as the whole clustered dataset. The whole annotated dataset can subsequently be used in downstream quantification and statistical analysis.


To demonstrate this strategy, we applied clustering and DR to a dataset of brain cells from mock- or WNVinfected mice (Figure 3A). The whole dataset was clustered using FlowSOM (Figure 3B), and the clustered data were proportionally subsampled for plotting with UMAP (Figure 3C). In this case, the number of cells extracted from each sample were in proportion to the total cells recovered from each brain, ensuring that DR plots accurately reflected sample composition (Figure 3D). By examining expression heatmaps (Figure 3E) and colored DR plots (Figure 3F) we manually determined cluster population identities, and these annotations were applied to both the subsampled and full datasets.

400 The choice of markers used to inform the generation of clusters/DR results is dependent on the overall goal 401 of analysis. Cellular markers may be broadly categorized into two groups: static (stably expressed) or 402 dynamic (changing expression). Typically, statically expressed markers are helpful for identifying consistent cell types (e.g. T cells, B cells, etc.), whereas dynamic markers are helpful for identifying reactive 403 404 cellular states (e.g. activated, resting, etc.). When seeking to discover novel cellular populations or states, 405 incorporating all cellular markers may be of benefit, as both stable cell types and dynamic cellular states 406 will be captured in separate clusters. However, a more selective approach may also be desired, such as using statically expressed markers to capture known populations within clusters, and then examining those stable 407 408 populations for dynamic changes in activation status. Along with domain-specific knowledge about the 409 expression patterns of various markers, it is possible to identify markers that contribute most to the level of 410 variance across a given dataset. To do this we can using principal component analysis (PCA), (Figure 3G-H) to determine the relative contribution of each cellular marker to the overall variance of the dataset. 411

412 Multi-level immune profiling

Figurice 3 int doi: https://doi.org/10.1101/2020.10.22.349563; this version posted October 23, 2020. The copyright holder for this preprint (which was not certified by peer review) is the author/funder, who has granted bioRxiv a license to display the preprint in perpetuity. It is made **A** CNS dataset **B** FlowSOM
available unger @MapBY-NC-ND 4.0 Interpatiented CNS
WNV-infected CNS








Figure 3. Clustering and dimensionality reduction using Spectre

A) A dataset of cells isolated from mock- or WNV-infected CNS were used to demonstrate clustering and DR in Spectre. B) FlowSOM clustering was performed on the full dataset, which was C) then subsampled and plotted using UMAP. D) Parsing the dataset by each experimental group reveals substantial changes to immune populations. E-F) An examination of marker expression on each cluster allows for a user-determined annotation into biological relevant cell types. G) Analysis using PCA allows for visualization of the data variance, and H) the relative contribution of markers to the first two principle components.

413 Many populations of interest, such as hematopoietic stem cells (HSC) in the BM, are of very low frequency 414 within individual samples. As such, representation of these rare subsets may be extremely sparse on DR plots, relative to more abundant populations. Moreover, more nuanced clustering and analyses of such 415 416 populations are often desired, but the global data structures provided by more abundant subsets of cells may 417 dominate the analysis. To address this, we extended our analysis approach to enable the exploration of data 418 at multiple levels. Upon clustering the complete dataset, clusters representing rare or novel populations can be isolated and re-clustered independent of other populations and thereafter annotated in greater detail. 419 Expanding on this approach, multiple lineages can be split and profiled independently, then re-merged, 420 421 retaining detailed cluster annotations for combined plotting and quantitative analysis. This process is 422 conceptually similar to hierarchical approaches such as hierarchical-SNE (h-SNE, [52]), but accommodates 423 a more bespoke tailoring of the analytical process.

To demonstrate this, we examined a dataset of BM HSCs generated by mass cytometry. Clusters containing HSC and progenitors (denoted by CD117 expression) (Figure 4A) were extracted from the full dataset, and subjected to independent clustering and DR (Figure 4B). This independent analysis allowed for a more detailed assessment of low frequency subsets (Figure 4C) that were not easily assessed in the plots from the full dataset.

429 Automated cellular classification and label transfer between aligned datasets

430 A crucial application of computational analysis is *discovery* – defining novel subsets and/or investigating 431 experimental changes in novel subsets or states in new diseases, tissues, or experimental conditions. Analytical approaches to this often depend heavily on unsupervised techniques, such as clustering and DR. 432 Such analyses culminate in an annotated datasets in which each cell is manually assigned a population label 433 434 by the user. Alternatively, other studies seek to apply a more *repetitive* analytical process, using semi-435 supervised tools, such as automatic gating, to replicate a method of population identification over a large number of samples. While effective, many of these tools rely on some form of gating strategy and forgo the 436 use of any unsupervised techniques to expedite gating, thereby limiting the possible identification of new 437




Figure 4. Multi-level analysis for profiling of rare populations

A) A UMAP plot (left) where clusters representing stem cell and progenitor subsets were identified via CD117 expression. Through cross-referencing against FlowSOM clusters (right), these cells were B) subjected to new clustering, subsampling, and plotting using UMAP. C) Expression color plots reveal low frequency cellular subsets that were difficult to otherwise detect on the full UMAP plot.

438 or complex/overlapping populations. In contrast, machine learning-based approaches provide the 439 opportunity for automated transfer of cellular labels from an annotated dataset to a novel dataset, following alignment of the datasets into a single feature space (as demonstrated in Figure 2). To facilitate this, we 440 provided functions within Spectre to train and run classifiers, a type of machine learning approach designed 441 442 to predict the class of given data points. As opposed to clustering, which groups cells together based on marker expression, classifiers assign cells a label based on 'training' data. This training data could be 443 previously gated, clustered, or annotated and is used by the classifier to determine how given features 444 (marker expression) relates to a class (cluster or cell phenotype). 445

Inspired in part by Seurat's mutual nearest neighbors approach for data integration [19, 20], we reasoned 446 447 that a simple nearest neighbors approach could facilitate rapid label transfer between datasets. In Spectre, 448 we implemented a k-Nearest Neighbor (kNN) classifier, which classifies unlabeled cells based on the label of their Kth nearest neighbor within the training dataset. To determine the accuracy of the kNN classifier, 449 450 we divided a labelled dataset into two halves (Supplementary Figure 7A) - one half retained cellular labels 451 and served as the 'training' dataset, and the other half had their cellular labels hidden, and served as the validation dataset. The kNN model was trained on the labelled data and applied to the unlabeled data 452 (Supplementary Figure 7B), to compare against their original cellular labels (Supplementary Figure 7C). 453 454 Generally, a kNN classification scheme with k = 1 tended to predict cellular labels between datasets with high accuracy (> 98.6%), although we found that increasing k here marginally increased this accuracy to 455 over 98.9% in this case (Supplementary Figure 7D). 456

To demonstrate the kNN classification process in an experimental context, we aligned two datasets of BM cells using CytoNorm (shown in Figure 2D-E), where the first dataset (batch A) contained labelled cells, and the second dataset (batch B) did not (Figure 5A). Following alignment, the data was split into the labeled and unlabeled datasets (Figure 5B), where the kNN classifier was trained on the labelled cells (batch A) and applied to unlabeled cells (batch B) with k = 1 (Figure 5C). When we compared the predicted cellular labels to manually annotated cellular labels, we found the classifier was able to accurately transfer cellular




Figure 5. Cellular classification and label transfer

A) A dataset of two batches of BM cells, following alignment with CytoNorm. One batch contains annotated clusters, and the other does not. B) The dataset is split into each batch (labeled and unlabeled), and a kNN classifier trained on the labelled dataset with k = 1. C) The kNN classifier was then applied to the unlabeled dataset, D) resulting in an accurate transfer of cellular labels between datasets.

labels between the two datasets (Figure 5D). Merging of the two datasets the resulted in a fully annotateddataset (Figure 5E).


465 Quantitative and statistical analysis


The endpoint of most analytical workflows is to make quantitative and statistical comparisons between 466 experimental groups. Many components in this workflow can be automated to simplify the processing steps 467 468 and reduce the time taken to generate relevant statistics and plots. To facilitate this, we have developed a 469 series of functions to summarize a dataset rapidly at either the cluster or cellular population level, resulting in a series of summary tables. For each population in each sample, these tables summarize the proportion 470 471 of cells, total cell counts, marker expression levels, and proportion of cells that are 'positive' for each 472 cellular marker, which can then be used to generate quantitative plots. The generation of grouped scatter or violin plots (Figure 6A) provide a simple method to assess changes of a single feature (e.g. number of 473 474 infiltrating macrophages per brain) between experimental groups, including grouped or pairwise statistical comparisons. However, more global statistical analyses are often desired. The generation of z-score 475 476 heatmaps (Figure 6B) provide an overview of relative changes between samples, with optional clustering on samples (rows) or features (columns) based on similarity. Additionally, the result of pairwise 477 478 comparisons between groups can be indicated for each heatmap column, revealing statistical significance 479 for uncorrected or false discovery rate (FDR)-corrected p-values. Furthermore, PCA plots (Figure 6C) and 480 volcano plots (Figure 6D) provide a further global view of how these immune features differentiate samples within the experimental context. 481

482

483 **DISCUSSION**

484 Challenges in cytometric analysis





Total = 48 variables

Figure 6. Quantification and statistical analysis

Quantitative and statistical analysis of CNS samples. A) Grouped dot plots with underlying violin plots indicating the relative cells per sample for each population in the dataset. Pairwise statistical analysis for non-parametric data was performed using a Mann-Whiteny/Wilcox test, where p < 0.05 was considered significant. Importantly, data for each plot were not adjusted for multiple comparisons. B) The z-score for each column of data was calculated, and plotted using the make.pheatmap function. Each row represents a sample, and each column represents a measured immune feature. Rows and columns were clustered using a uclidean distance. Pairwise comparisons between experimental groups were calculated using a Mann Whiteny/Wilcox test, and the significance results plotted as significant (p < 0.05) or non-significant (p > 0.05) for each column. To adjust for multiple comparisons, p-value results were corrected using a FDR correction, and results plotted as significant (p < 0.05) or non-significant (p > 0.05) for each column. C) PCA results show the distribution of each sample, and D) volcano plots show the relative fold-change increase or decrease of each immune feature (X axis), and the inverse p-value of each change (Y axis).

485 Historically, analysis of cytometry data involves a manual process of identifying cells, usually through 486 drawing a gate on cells plotted with one marker vs another. As the number of features that can be measured on cells has continued to increase, the process of gating alone has become untenable for comprehensive 487 488 analysis of large and complex datasets. This is in part due to the sheer number of gates that need to be drawn 489 to fully parse a dataset, and in part due to the bias that is introduced by selectively including/excluding cells 490 using one or two markers at a time sequentially. As a result, a number of computationally-driven approaches to data analysis have been developed, including automated gating, clustering, DR, and classification. These 491 enable automated, data-driven data processing that may be performed in an unsupervised or semi-492 493 supervised fashion. However, these approaches come with their own limitations, including operational 494 constraints, such as the number of cells that can be analyzed, given the available computational resources, 495 the possibility of inaccurate identification of specific biologically-relevant populations in the dataset, and a 496 lack of reproducibility between analytical runs (particularly when applying to new datasets). Critically, 497 many analysis tools are developed as standalone packages that operate on specific data formats. This limits 498 the interoperability of tools, especially when attempting to combine tools developed in different fields (such as cytometry and scRNAseq). As a solution, we developed Spectre: an adaptable and easy-to-use package 499 for analyzing high-dimensional cytometry data. Spectre enhances existing computational tools through 500 501 strategic implementation and customization of high-performance computational tools, and provision of 502 wrapper functions to simplify and improve their flexibility and interoperability. This package expanded on 503 many of the key aspects of the CAPX workflow [32] that has been utilized in a number of studies [53-57].

504 As the foundation for data in Spectre, data.table allows for easy handling of large datasets

505 Many cytometry or single-cell analysis tools operate on custom data formats, such as the flowFrame [35], 506 SingleCellExperiment [21], or Seurat objects [19]. Each custom data format may contain numerous 507 elements, made up of both base R and custom data formats. Typically, primary cellular expression data is 508 stored in a table or matrix, and separate elements contain metadata for each cell (cell number, sample, 509 group, etc.,) or feature (marker names, parameters names, voltages, gene numbers, etc). Importantly, 510 manipulated data (such as ArcSinh transformations or normalization) may be contained within the primary 511 data table/matrix, or structured separately. The complexity of these custom formats, and the significant 512 differences in structure between them, results in difficulties in converting one data format to another. While 513 custom data formats may be convertible into other formats, including conversion to one of R's base data 514 formats, these conversions are non-trivial, and often result in the loss of important metadata. Of note, 515 cytometry data is typically structured with rows representing cells, and columns representing cellular 516 features, though this is transposed in scRNAseq data.

517 The foundation of Spectre is built upon the *data.table* package, which enhances R's base *data.frame* format. The efficiency of *data.table* in performing basic data manipulation operations, such as filtering, ordering, 518 519 importing, and exporting, makes it suitable for processing large HD cytometry datasets. While there are R 520 packages which operate on *data.frame* formats (e.g. dplyr [58]), our simple benchmarking measurements show *data.table* to be faster when handling large datasets (Supp. Figure 1). The simple tabular format allows 521 522 for interoperability between Spectre, basic R functions, and other functions from cytometry or single-cell-523 specific packages, by storing all the relevant information for each cell in a single, high-performance table. Here, columns desired for use with each function can be easily specified, and new columns that are added 524 as a result of analysis are easily identifiable through the use of regular expression patterns (such as 525 526 "CD4 asinh", "CD4 aligned", etc.). Column metadata in this context is less relevant for processing and 527 analysis, but can still be imported and managed with the main data should the user choose to do so. By 528 designing all our functions to operate on this simple *data.table* structure, we remove the requirements for 529 users to convert their data into specific formats for different functions, greatly improving usability.

530 Usage flexibility

We have demonstrated the use of Spectre through various workflows. While they each illustrate an approach to analysis, they are flexible in such a way that functions can be run in any order, and can be replaced with others (including those that are not built into Spectre). For instance, it is possible to run a DR tool on the dataset before clustering, or replacing Spectre's run kNN classifier function with another

535 function (like a Decision Tree classifier [59]). Additionally, most of Spectre's functions are designed such 536 that they do not rely on each other to operate. A function may simply be executed on its own without having 537 to run others beforehand. Some functions require data to be in a particular layout with appropriate calculations (e.g. calculating MFI across samples for heatmap creation with make.pheatmap), but this can 538 539 be achieved either with (using write.sumtables) or without (using other functions or manual creation) 540 Spectre. In addition, Spectre's functions are equipped with parameters that allow users to customize how they operate. The majority of these parameters come with a set of default values that are based on the 541 original implementation of the tools. The capacity to run Spectre's functions with either default or custom 542 543 parameter values make Spectre both simple and customizable.

544 Docker for accessibility

Spectre provides wrapper functions to many existing computational tools, and thus require users to install 545 546 the corresponding R packages before using them. Inadvertently, this creates a cycle of package dependencies as these packages often rely on other existing R packages. Managing such dependencies has 547 548 proven to be a challenge even to advanced analysts, as there exists a myriad of ways to install, remove, and 549 update packages (e.g. CRAN, Bioconductor [60], devtools [61]). Moreover, some packages rely on users 550 to have other softwares (e.g. Xcode in Mac, Rtools in Windows) or compilers available on their computer. 551 This may pose as a major hurdle for those who are not familiar with R. To address this, we made Spectre 552 available as a Docker image. The image is a prepared environment loaded with RStudio for users to interact with R code (write and run) and all the libraries required by Spectre. By downloading the Docker image 553 554 and launching it as a self-contained computing environment, users will be able to run their analysis without going through the complicated setup process. While Docker introduces an additional layer between user's 555 556 physical computing resources and the analytical tools, previous work by IBM indicates the performance degradation to be negligible [62]. 557

558 Versatility of application

559 The structure of all forms of cytometry data, including flow, spectral or mass, are inherently similar: each row is an individual cell, whilst each column is an individual marker or feature measured on those cells. 560 561 The values in such a table indicate the signal intensity for each individual marker within the panel on each cell. Thus, Spectre can be used on datasets generated by both flow (including spectral) and mass cytometry, 562 563 following compensation (or spectral unmixing), and initial cleanup gating. As a result, Spectre can also be used on other forms of cellular data, such single-cell RNA seq data, following some additional pre-564 processing steps. Additionally, Spectre can be used to analyze HD imaging data, such as that generated by 565 IMC, once cellular segmentation has been performed. 566

567

568 DATA AVAILABILITY

569 Spectre code and demonstration data used in this paper are available at

- 570 <u>https://github.com/immunedynamics/spectre</u>.
- 571

572 ACKNOWLEDGEMENTS

This work was supported by a grant from the Marie Bashir Institute for Infectious Disease and Biosecurity. 573 This work was supported in part by NH&MRC Project grant 1088242 and grant from the Merridew 574 575 Foundation to N.J.C.K. Authors T.M.A. and F.M-W. are supported by the International Society for the 576 Advancement of Cytometry (ISAC) Marylou Ingram Scholars program. G.H.P. is supported by the Australian Government Research Training Program (RTP) Scholarship. A.G.S. is supported by the 577 Australian Government RTP Scholarship and The University of Sydney Postgraduate Merit Award. We 578 would like to acknowledge the contribution of our colleagues who inspired the design and implementation 579 of Spectre. 580

581 **REFERENCES**

- Bendall, S.C., et al., Single-cell trajectory detection uncovers progression and regulatory
 coordination in human B cell development. Cell, 2014. 157(3): p. 714-25.
- 584 2. Mair, F., et al., *The end of gating? An introduction to automated analysis of high dimensional*585 *cvtometry data.* Eur J Immunol, 2016. 46(1): p. 34-43.
- Brinkman, R.R., *Improving the Rigor and Reproducibility of Flow Cytometry-Based Clinical Research and Trials Through Automated Data Analysis.* Cytometry A, 2020. 97(2): p. 107-112.
- Levine, J.H., et al., *Data-Driven Phenotypic Dissection of AML Reveals Progenitor-like Cells that Correlate with Prognosis.* Cell, 2015. 162(1): p. 184-97.
- 5. Van Gassen, S., et al., *FlowSOM: Using self-organizing maps for visualization and interpretation of cytometry data.* Cytometry A, 2015. 87(7): p. 636-45.
- 592 6. Samusik, N., et al., *Automated mapping of phenotype space with single-cell data*. Nat Methods,
 593 2016. 13(6): p. 493-6.
- van der Maaten, L. and G. Hinton, *Visualizing Data using t-SNE*. Journal of Machine Learning
 Research, 2008. 9: p. 2579-2605.
- van der Maaten, L., *Accelerating t-SNE using Tree-Based Algorithms*. Journal of Machine Learning
 Research, 2014. 15: p. 3221-3245.
- McInnes, L., J. Healy, and J. James Melville, UMAP: Uniform Manifold Approximation and
 Projection for Dimension Reduction. arXiv, 2018: p. 1802.03426.
- Setty, M., et al., *Wishbone identifies bifurcating developmental trajectories from single-cell data*.
 Nat Biotechnol, 2016. 34(6): p. 637-45.
- Li, H., et al., *Gating mass cytometry data by deep learning*. Bioinformatics, 2017. 33(21): p. 34233430.
- 604 12. Chen, Y., et al., *Single-cell classification using learned cell phenotypes*. bioRxiv, 2020: p.
 605 2020.07.22.216002.
- Kaushik, A., et al., *CyAnno: A semi-automated approach for cell type annotation of mass cytometry datasets.* bioRxiv, 2020: p. 2020.08.28.272559.
- 608 14. Chen, H., et al., *Cytofkit: A Bioconductor Package for an Integrated Mass Cytometry Data Analysis*609 *Pipeline*. PLoS Comput Biol, 2016. **12**(9): p. e1005112.
- Bruggner, R.V., et al., *Automated identification of stratifying signatures in cellular subpopulations*.
 Proc Natl Acad Sci U S A, 2014. 111(26): p. E2770-7.
- 612 16. Chevrier, S., et al., *Compensation of Signal Spillover in Suspension and Imaging Mass Cytometry*.
 613 Cell Syst, 2018. 6(5): p. 612-620 e5.

- 614 17. Nowicka, M., et al., *CyTOF workflow: differential discovery in high-throughput high-dimensional*615 *cytometry datasets.* F1000Res, 2017. 6: p. 748.
- 616 18. Weber, L.M., et al., *diffcyt: Differential discovery in high-dimensional cytometry via high-*617 *resolution clustering.* Commun Biol, 2019. 2: p. 183.
- Butler, A., et al., Integrating single-cell transcriptomic data across different conditions,
 technologies, and species. Nat Biotechnol, 2018. 36(5): p. 411-420.
- 620 20. Stuart, T., et al., *Comprehensive Integration of Single-Cell Data*. Cell, 2019. 177(7): p. 1888-1902
 621 e21.
- 622 21. Amezquita, R.A., et al., *Orchestrating single-cell analysis with Bioconductor*. Nat Methods, 2020.
 623 17(2): p. 137-145.
- 22. Zeng, C., et al., *Pseudotemporal Ordering of Single Cells Reveals Metabolic Control of Postnatal beta Cell Proliferation*. Cell Metab, 2017. 25(5): p. 1160-1175 e11.
- 626 23. Tran, H.T.N., et al., A benchmark of batch-effect correction methods for single-cell RNA
 627 sequencing data. Genome Biol, 2020. 21(1): p. 12.
- 628 24. Regev, A., et al., *The Human Cell Atlas*. Elife, 2017. 6.
- 629 25. Korsunsky, I., et al., *Fast, sensitive and accurate integration of single-cell data with Harmony*. Nat
 630 Methods, 2019. 16(12): p. 1289-1296.
- 631 26. Hao, Y., et al., *Integrated analysis of multimodal single-cell data*. bioRxiv, 2020: p.
 632 2020.10.12.335331.
- Lin, Y., et al., scMerge leverages factor analysis, stable expression, and pseudoreplication to
 merge multiple single-cell RNA-seq datasets. Proc Natl Acad Sci U S A, 2019. 116(20): p. 97759784.
- Bahne, F., et al., *Per-channel basis normalization methods for flow cytometry data*. Cytometry A,
 2010. 77(2): p. 121-31.
- 638 29. Schuyler, R.P., et al., *Minimizing Batch Effects in Mass Cytometry Data*. Front Immunol, 2019. 10:
 639 p. 2367.
- 640 30. Van Gassen, S., et al., *CytoNorm: A Normalization Algorithm for Cytometry Data*. Cytometry A,
 641 2020. 97(3): p. 268-278.
- 642 31. Trussart, M., et al., *Removing unwanted variation with CytofRUV to integrate multiple CyTOF*643 *datasets.* Elife, 2020. 9.
- 644 32. Ashhurst, T.M., et al., *Analysis of the Murine Bone Marrow Hematopoietic System Using Mass and*645 *Flow Cytometry*. Methods Mol Biol, 2019. **1989**: p. 159-192.
- 646 33. Niewold, P., et al., *Evaluating spectral cytometry for immune profiling in viral disease*. Cytometry
 647 A, 2020.

- 648 34. Dowle, M. and A. Srinivasan, data.table: Extension of `data.frame`. R package version 1.13.0.
 649 https://CRAN.R-project.org/package=data.table. 2020.
- 35. Hahne, F., et al., *flowCore: a Bioconductor package for high throughput flow cytometry.* BMC
 Bioinformatics, 2009. 10: p. 106.
- 652 36. Van Gassen, S., B. Callebaut, and Y. Saeys, *FlowSOM: Using self-organizing maps for*653 *visualization and interpretation of cytometry data.*654 *http://bioconductor.org/packages/release/bioc/html/FlowSOM.html.* 2020.
- 655 37. Krijthe, J.H., *Rtsne: T-Distributed Stochastic Neighbor Embedding using a Barnes-Hut*656 *Implementation. <u>https://github.com/jkrijthe/Rtsne</u>. 2015.*
- 657 38. Konopka, T., *umap: Uniform Manifold Approximation and Projection. R package version 0.2.5.0.* 658 <u>https://CRAN.R-project.org/package=umap</u>. 2020.
- 659 39. Team, R.C., *R: A Language and Environment for Statistical Computing*. <u>https://www.R-</u>
 660 project.org/. 2020.
- 40. Wickham, H., ggplot2: Elegant Graphics for Data Analysis. <u>https://ggplot2.tidyverse.org</u>. 2016:
 Springer-Verlag New York.
- 41. Kolde, R., pheatmap: Pretty Heatmaps. R package version 1.0.12. <u>https://CRAN.R-</u>
 664 project.org/package=pheatmap. 2019.
- 42. Kassambara, A., ggpubr: 'ggplot2' Based Publication Ready Plots. R package version 0.4.0.
 https://CRAN.R-project.org/package=ggpubr. 2020.
- 667 43. Blighe, K., S. Rana, and M. Lewis, *EnhancedVolcano: Publication-ready volcano plots with*668 *enhanced colouring and labeling. R package version 1.6.0.*669 *https://github.com/kevinblighe/EnhancedVolcano.* 2020.
- 670 44. Kuhn, M., caret: Classification and Regression Training. R package version 6.0-86.
 671 https://CRAN.R-project.org/package=caret. 2020.
- 672 45. Beygelzimer, A., et al., FNN: Fast Nearest Neighbor Search Algorithms and Applications. R
 673 package version 1.1.3. <u>https://CRAN.R-project.org/package=FNN</u>. 2019.
- 46. Parks, D.R., M. Roederer, and W.A. Moore, *A new "Logicle" display method avoids deceptive effects of logarithmic scaling for low signals and compensated data*. Cytometry A, 2006. 69(6): p.
 541-51.
- 47. Moore, W.A. and D.R. Parks, Update for the logicle data scale including operational code
 678 implementations. Cytometry A, 2012. 81(4): p. 273-7.
- 679 48. Bendall, S.C., et al., *Single-cell mass cytometry of differential immune and drug responses across*680 *a human hematopoietic continuum*. Science, 2011. 332(6030): p. 687-96.

- 49. Weber, L.M. and M.D. Robinson, *Comparison of clustering methods for high-dimensional single-*
- 682 *cell flow and mass cytometry data.* Cytometry A, 2016. **89**(12): p. 1084-1096.
- 50. Linderman, G.C., et al., *Fast interpolation-based t-SNE for improved visualization of single-cell RNA-seq data.* Nat Methods, 2019. 16(3): p. 243-245.
- 51. Belkina, A.C., et al., Automated optimized parameters for T-distributed stochastic neighbor
 embedding improve visualization and analysis of large datasets. Nat Commun, 2019. 10(1): p.
 5415.
- van Unen, V., et al., *Visual analysis of mass cytometry data by hierarchical stochastic neighbour embedding reveals rare cell types.* Nat Commun, 2017. 8(1): p. 1740.
- 53. Vuckovic, S., et al., *Inverse relationship between oligoclonal expanded CD69- TTE and CD69+ TTE cells in bone marrow of multiple myeloma patients.* Blood Adv, 2020. 4(19): p. 4593-4604.
- 692 54. Marsh-Wakefield, F., et al., *IgG3 (+) B cells are associated with the development of multiple*693 *sclerosis.* Clin Transl Immunology, 2020. 9(5): p. e01133.
- 694 55. Marsh-Wakefield, F., et al., *Mass Cytometry Discovers Two Discrete Subsets of CD39(-)Treg*695 Which Discriminate MGUS From Multiple Myeloma. Front Immunol, 2019. 10: p. 1596.
- 56. Shinko, D., et al., *Mass Cytometry Reveals a Sustained Reduction in CD16(+) Natural Killer Cells Following Chemotherapy in Colorectal Cancer Patients.* Front Immunol, 2019. 10: p. 2584.
- 698 57. Hayashida, E., et al., Zika virus encephalitis in immunocompetent mice is dominated by innate
 699 immune cells and does not require T or B cells. J Neuroinflammation, 2019. 16(1): p. 177.
- 58. Wickham, H., et al., *dplyr: A Grammar of Data Manipulation. R package version 0.8.5.* https://CRAN.R-project.org/package=dplyr. 2020.
- 59. Belson, W.A., "*Matching and prediction on the principle of biological classification*." Journal of
 the Royal Statistical Society: Series C (Applied Statistics) 8.2, 1959: p. 65-75.
- Morgan, M., BiocManager: Access the Bioconductor Project Package Repository. R package
 version 1.30.10. https://CRAN.R-project.org/package=BiocManager. 2019.
- Wickham, H., J. Hester, and W. Chang, *devtools: Tools to Make Developing R Packages Easier. R package version 2.3.0. <u>https://CRAN.R-project.org/package=devtools</u>. 2020.*
- Felter, W., et al., *An updated performance comparison of virtual machines and Linux containers*,.
 2015 IEEE International Symposium on Performance Analysis of Systems and Software (ISPASS),
 2015: p. 171-172.
- 711

712

713 FIGURE LEGENDS

714 Figure 1. Spectre analysis overview

715 An overview of Spectre's analysis workflow. A) Data preparation steps including sample, group, and batch 716 annotation, in addition to ArcSinh transformation. B) Batch alignment using CytoNorm. C) Clustering and 717 dimensionality reduction, along with marker expression plotting and expression heatmaps. D) 718 Ouantification and statistical analysis through z-score/fold-change heatmaps and grouped dot plots. E) 719 Application of Spectre's analysis workflow to data generated by different technologies, including analyzing 720 a split murine spleen sample by flow cytometry, spectral cytometry, and mass cytometry. Also shown are 721 unrelated PBMC data analyzed by single-cell RNA sequencing (scRNAseq) acquired from 10X genomics, 722 via the Seurat webpage (https://satijalab.org/seurat/v3.2/pbmc3k_tutorial.html) and imaging data generated 723 by Imaging Mass Cytometry (IMC).

724 Figure 2. Batch alignment using CytoNorm

725 Batch alignment process using CytoNorm. A) Reference samples acquired with each batch are extracted 726 from the *data.table*, and B) clustered using FlowSOM. Metacluster-specific quantile conversion models are 727 then calculated. C) Cells from all samples/batches are mapped to the FlowSOM grid, and assigned to their nearest metacluster. Cells are then aligned using the metacluster-specific quantile conversion models D) 728 729 Two sets of BM samples with synthetic batch effects introduced in a population-specific manner. E) 730 CytoNorm alignment was performed using a metacluster-specific alignment process (fine alignment). F) 731 Two sets of BM samples generated with slightly different panels, but targeting the same cellular markers, resulting in significant batch effects. G) CytoNorm is initially performed on the whole dataset (coarse 732 733 alignment) by mapping the entire dataset into a single metacluster, where H) subsequent FlowSOM 734 clustering allowed for further metacluster-specific alignment (fine alignment).

735 Figure 3. Clustering and dimensionality reduction using Spectre

A) A dataset of cells isolated from mock- or WNV-infected CNS were used to demonstrate clustering and
DR in Spectre. B) FlowSOM clustering was performed on the full dataset, which was C) then subsampled
and plotted using UMAP. D) Parsing the dataset by each experimental group reveals substantial changes to
immune populations. E-F) An examination of marker expression on each cluster allows for a userdetermined annotation into biological relevant cell types. G) Analysis using PCA allows for visualization
of the data variance, and H) the relative contribution of markers to the first two principle components.

742 Figure 4. Multi-level analysis for profiling of rare populations

A) A UMAP plot (left) where clusters representing stem cell and progenitor subsets were identified via
CD117 expression. Through cross-referencing against FlowSOM clusters (right), these cells were B)
subjected to new clustering, subsampling, and plotting using UMAP. C) Expression color plots reveal low
frequency cellular subsets that were difficult to otherwise detect on the full UMAP plot.

747 Figure 5. Cellular classification and label transfer

A) A dataset of two batches of BM cells, following alignment with CytoNorm. One batch contains annotated clusters, and the other does not. B) The dataset is split into each batch (labeled and unlabeled), and a kNN classifier trained on the labelled dataset with k = 1. C) The kNN classifier was then applied to the unlabeled dataset, D) resulting in an accurate transfer of cellular labels between datasets.

752 Figure 6. Quantification and statistical analysis

Quantitative and statistical analysis of CNS samples. A) Grouped dot plots with underlying violin plots indicating the relative cells per sample for each population in the dataset. Pairwise statistical analysis for non-parametric data was performed using a Mann-Whiteny/Wilcox test, where p < 0.05 was considered significant. Importantly, data for each plot were not adjusted for multiple comparisons. B) The z-score for each column of data was calculated, and plotted using the make.pheatmap function. Each row represents a sample, and each column represents a measured immune feature. Rows and columns were clustered using euclidean distance. Pairwise comparisons between experimental groups were calculated using a Mann

760	Whiteny/Wilcox test, and the significance results plotted as significant (p < 0.05) or non-significant (p $>$
761	0.05) for each column. To adjust for multiple comparisons, p-value results were corrected using a FDR
762	correction, and results plotted as significant ($p < 0.05$) or non-significant ($p > 0.05$) for each column. C)
763	PCA results show the distribution of each sample, and D) volcano plots show the relative fold-change
764	increase or decrease of each immune feature (X axis), and the inverse p-value of each change (Y axis).
765	
766	
767	
768	
769	
770	
771	
772	
773	
774	
775	
776	
777	
778	
779	
780	

781 Supplementary Figure 1. Data processing speeds

- A) Benchmarking of files reading/writing with data.table or base R functions. B) Benchmarking of common
- data operations using data.table or Dplyr. 12,108,595 cells with approx. 25 dimensions were analyzed on a
- virtual windows machine with 16 GB of RAM.

785 Supplementary Figure 2. Data visualization and ArcSinh transformation

A) Brief comparison of plot visualization options as implemented in FlowJo, with the axis set as linear
(with default maximum values), linear (with adjusted maximums), logarithmic, or bi-exponential/logicle.
B) Direct comparison of raw data using a bi-exponential/logicle axis with data subject to ArcSinh
transformation and plotted on a linear axis.

790 Supplementary Figure 3. ArcSinh transformation of mass cytometry data

A) Plots of raw data in FlowJo with various width basis (WB) settings for the x-axis (168Er CD8a).
Minimum x-axis value for each plot was adjusted to the highest possible value, within the restrictions
imposed by FlowJo. Y-axis width basis was fixed at -100. Both x- and y-axis positive decades were set at
4.75. B) Plots of ArcSinh transformed data from R with various co-factors used for the marker on the xaxis. Y-axis co-factor was fixed at 10.

796 Supplementary Figure 4. ArcSinh transformation of conventional and spectral cytometry data

Plots of data generated by a conventional (A-B) and spectral (C-D) flow cytometry system. The conventional flow cytometry system was a 10-laser BD LSR-II, and the spectral flow cytometry system was a 5-laser Cytek Aurora. A) Plots with of conventional flow cytometry raw data in FlowJo with various width basis (WB) settings for the x-axis (BUV805 CD8a). Minimum x-axis value for each plot was adjusted to the highest possible value, within the restrictions imposed by FlowJo. Y-axis width basis was fixed at – 100. Both x- and y-axis positive decades were set at 4.29. B) Plots of ArcSinh transformed data from R with various co-factors used for the marker on the x-axis. Y-axis co-factor was fixed at 1000. C) Plots of

spectral cytometry raw data in FlowJo with various width basis (WB) settings for the x-axis (BUV805 CD8a). Minimum x-axis value for each plot was adjusted to the highest possible value, within the restrictions imposed by FlowJo. Y-axis width basis was fixed at -631. Both x- and y-axis positive decades were set at 5.68. D) Plots of ArcSinh transformed data from R with various co-factors used for the marker on the x-axis. Y-axis co-factor was fixed at 5000.

809 Supplementary Figure 5. ArcSinh, normalized, and noise reduced data

810 Transformations and manipulations performed on conventional flow cytometry data. A) Raw data plot of

811 FITC Ly6C vs BUV805 CD8a. B) ArcSinh transformed data (a-axis) using a co-factor (CF) of 1000. C)

812 ArcSinh transformed data (x-axis) redistributed between 0 and 1 using the *do.normalise* function. D)

813 ArcSinh transformed data (a-axis) with data points below 0 changed to 0 using the *do.noise.reduce* function,

to reduce the distribution of negative cells below zero.

815 Supplementary Figure 6. Dimensionality reduction crowding and computation time

A-B) Plots of tSNE analysis performed with different numbers of total cells, colored by relative density, or

817 cell type. C) Computation time for tSNE on different sized datasets, plotted on a logarithmic scale.

818 Supplementary Figure 7. Evaluating kNN classification accuracy

A) A single labelled dataset was split into halves – one half retaining population labels (training data) and one half having population labels hidden (validation data). B) A kNN model was trained on the training data and applied to the validation data. C) Comparison of the original labels on the training data and the predicted labels on the validation data. D) A graph showing the accuracy for the kNN classifier for different values of *k*.